

Extraction Forceps Function

Dental forceps are pliers-like instruments that are used to grasp teeth.

A dentist will usually have a number of different ones on hand, each having a design that's tailored to:
The specific shape of the tooth it's intended to remove (such as incisors vs. molars having two roots vs. molars having three roots).
The location of the tooth in the mouth (such as front vs. back, or left lower molar vs. right upper molar.) How they're used.

A dentist will grasp a tooth with their forceps and then firmly and deliberately rock it back and forth as much as it will.

Because the bone that surrounds the tooth is compressible, the socket will gradually start to expand.

As it does, the range of the dentist's side-to-side motions will increase.

In addition to this rocking motion, the dentist will also rotate the tooth back and forth.

This twisting action helps to rip and tear the tooth away from the ligament that binds it in place.

At some point, the socket will be enlarged enough, and the ligament torn away enough, that the tooth can be easily removed.

Handle Selection For Extracting Forceps

Dear Valued Buyers these are or Handle Styles



English Pattern



English Pattern



English Pattern,
Anatomical Handle



American Pattern



Deep Grip Forceps
With Hole Pattern




Hole Pattern




Dear Valued Buyers if any Universal Fig. No. of Extracting Forceps not Found in our Catalogue you Can Refer the Fig. No. From These Diagram.
Dear Valued Buyers these are or Handle Styles

HERE IS THE UNIVERSAL FIG. NUMBERS IN BOX



HERE IS THE UNIVERSAL FIG. NUMBERS IN BOX



Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك



Article# IC6-001

Fig. 1
Upper Centrals
and Canines



Article# IC6-002

Fig. 2
Upper Laterals
and Canines



Article# IC6-003

Fig. 4
Lower Incisors
and Canines



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for
upper incisors and canines



Article# IC6-004

Fig. 7
Upper Premolars



Article# IC6-005

Fig. 8
Lower Premolars

■ Инструменты для реставрации
 ■ Instruments pour plombage
 ■ Füllungsinstrumente (mit Kontrawinkel)

■ Инструменти для реставрації
 ■ Instrumentos para obturaciones
 ■ ملء الصكوك



Article# IC6-006
Fig. 13
 Lower Premolars



Straight extracting forceps for upper incisors and canines



Article# IC6-007
Fig. 17
 Upper Molars, Right



Straight extracting forceps for upper incisors and canines



Article# IC6-008
Fig. 18
 Upper Molars, Left



Straight extracting forceps for upper incisors and canines



Article# IC6-009
Fig # 18A
 upper molars



Article# IC6-010
Fig # 19
 upper third molars



Article# IC6-011
Fig # 20
 lower third molar

Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك



Article# IC6-012
Fig. 21
Lower Molar



Article# IC6-013
Fig # 22
Lower Molar



Angled extracting forceps for upper roots



Article# IC6-014
Fig # 24
lower molar, right



Article# IC6-015
Fig # 29
upper roots



Article# IC6-016
Fig # 30
Upper Roots

■ Инструменты для реставрации
 ■ Instruments pour plombage
 ■ Füllungsinstrumente (mit Kontrawinkel)

■ Инструменти для реставрації
 ■ Instrumentos para obturaciones
 ■ ملء الصكوك

WISDOM EXTRACTING FORCEPS



Article# IC6-017
Routerier
(# 22 1/2 L)
lower molars and
wisdoms, left




Article# IC6-018
Routerier
(# 22 1/2 R)
lower molars and
wisdoms, right



Article# IC6-019
Routerier
(# 33 1/2 L)
lower molar and
wisdom roots, left



Article# IC6-020
Routerier
(# 33 1/2 R)
lower molar and
wisdom roots, right



Article# IC6-021
Fig. 186L
lower molars



Article# IC6-022
Fig. 186R
lower molars

Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك



Article# IC6-023

Fig # 31

Lower Roots Both Sides EF, Serrated, EU



Article# IC6-024

Fig # 32

Lower Roots Both Side EF, Serrated, EU



Article# IC6-025

33

lower roots



Article# IC6-026

33A

Lower Roots



Article# IC6-027

33M

Lower roots
with long beak

Инструменты для реставрации
Instruments pour plombage
Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації
Instrumentos para obturaciones
ملء الصكوك



Thin Tip

Article# IC6-028
Fig # 34
Upper incisors and
premolars
deep-gripping



Wide Tip

Article# IC6-029
Fig # 34A
Upper centrals / cains
deep-gripping



Thin Tip

Article# IC6-030
Fig # 35
upper premolars
deep-gripping



Wide Tip

Article# IC6-031
Fig # 35A
Upper incisors and
premolars
deep-gripping



Thin Tip

Article# IC6-032
Fig # 36
Lower premolars
deep-gripping



Wide Tip

Article# IC6-033
Fig # 36A
Lower premolars
deep-gripping



Article# IC6-034
Fig # F-7C
lower molar



Article# IC6-035
Fig # 336
lower premolars
deep-gripping

Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك



Article# IC6-036
Fig # 41
Upper Roots Both
Sides EF, Serrated, EU



Article# IC6-037
Fig # 44Z



Article# IC6-038
Fig # 46L
Lower Roots EF, Serrated EU



Article# IC6-039
Fig # 46XL
Apex Forceps Lower Root
Deep Root Pick



Article# IC6-040
Fig # 49
Upper Roots EF, Serrated, EU



Article# IC6-041
Fig # 51
Upper Roots



Article# IC6-042
Fig # 51A
Upper Roots



Article# IC6-043
Fig # 51AL



Article# IC6-044
Fig # 331
upper roots
deep-gripping

Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك



Article# IC6-045
Fig # 67
Upper Third Molars




Article# IC6-046
Fig # 67A
Upper Third Molars



Article# IC6-047
Fig # 67N
Lower third molars
deep-gripping



Article# IC6-048
Fig # 68
Lower roots



Article# IC6-049
Fig # 73
Lower Molars EF, Serrated, EU




Article# IC6-050
Fig # 73 A
lower molars



Article# IC6-051
Fig # 73 S



Article# IC6-052
Fig # 74
Lower Roots EF, Serrated, EU



Article# IC6-053
Fig # 74N
Lower Roots EF, Serrated, EU



Article# IC6-054
Fig # 75
Lower Premolars



Article# IC6-055
Fig # 79
Lower Third Molars

Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك



Article# IC6-056
Fig # 86A
Lower Molars



Article# IC6-057
Fig # 86B
lower molars



Article# IC6-058
Fig # 86C
lower molars



Article# IC6-059
Fig # 87
Lower Molars



Article# IC6-060
Fig # 89
Upper Molars,
Right



Article# IC6-061
Fig # 90
Upper Molars Left EF, Serrated, EU



Article# IC6-062
Fig # 114
Lower roots
and centrals



Article# IC6-063
Fig # 143
Lower roots
and centrals



Article# IC6-064
Fig # 144 (Balding)
Lower roots
and centrals



■ Инструменты для реставрации
 ■ Instruments pour plombage
 ■ Füllungsinstrumente (mit Kontrawinkel)

■ Інструменти для реставрації
 ■ Instrumentos para obturaciones
 ■ ملء الصكوك



Article# IC6-065
FIG. 97

Radici superiori molto sottili
Very thin upper roots
Racines supérieurs très fines
Raíces superiores muy delgada




Article# IC6-066
FIG. 151

Radici superiori
Upper roots
Racines supérieurs
Raíces superiores



Article# IC6-67
FIG. 351S

Radici superiori molto sottili
Very thin upper roots
Racines supérieurs très fines
Raíces superiores muy delgadas



Article# IC6-071
fig# 233
Lower Roots
EF, Serrated, EU



Article# IC6-068
fig # 97
Upper Roots EF,
Serrated, EU



Article# IC6-069
Fig. 130
(Ladmore, 190mm)
Upper wisdoms



Article# IC6-070
FIG. 333S

Radici inferiori molto sottili
Very thin lower roots
Racines inférieurs très fines



Article# IC6-072
146
Apex Forceps Lower Root
Deep Root Pick



Article# IC6-073
153
Apex Forceps Upper Root
Deep Root Pick



Article# IC6-074
149
Apex Forceps Upper Root
Deep Root Pick

Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك



Article# IC6-076

FIG. MEAD

Incisivi, canini, premolari e radici superiori
Upper incisors, canines, premolars and roots
Incisives, canines, prémolaires et racines supérieures
Incisivos, caninos, premolares y raíces superiores



Article# IC6-077

FIG. MEAD

1° - 2° - 3° molare - universale
1st - 2nd - 3rd molar - universal
1° - 2° - 3° molaire - universel
1° - 2° - 3° molare - universale



Article# IC6-078

FIG. MEAD

Incisivi e radici inferiori
Lower incisors and roots
Incisives et racines inférieures
Incisivos y raíces inferiores



Article# IC6-079




FIG. MEAD

1° - 2° - 3° molare inferiori
Lower 1st - 2nd - 3rd molar
1° - 2° - 3° molaire inférieurs



Article# IC6-075

Lower Molar
Fig. 23 Cow Horn



Article# IC6-080

Meissner
Extracting Forceps



Article# IC6-081

Meissner
Extracting Forceps



Article# IC6-082

Meissner
Extracting Forceps



Article# IC6-083

Comellas Forceps
Anatomical Handle

■ Инструменты для реставрации
 ■ Instruments pour plombage
 ■ Füllungsinstrumente (mit Kontrawinkel)

■ Инструменти для реставрації
 ■ Instrumentos para obturaciones
 ■ ملء الصكوك

(Deep Gripe Extracting Forceps)

A-TRAUMATIC EXTRACTION FORCEPS

Improved tapered beaks for complete root contact.
Cross serrations for a never slip grasp.
Easily enter the sub gingival.

- No Fracture
- No Tissue Damage
- Extraction With Easy



Article# IC6-085
Fig. 34A
Upper Anterior



Article# IC6-086
Fig. 35
Upper Universal



Article# IC6-087
Fig. 35A
Upper Universal



Article# IC6-088
Fig. 321
Lower Universal

Deep Root Pick Forceps



Article# IC6-089
Fig. 321A
Lower Root Molars



Article# IC6-090
Fig. 36
Lower Premolars
Deep-Gripping



Article# IC6-091
Fig. 36A
Lower Premolars
Deep-Gripping



Article# IC6-092
67A



Straight extracting forceps for upper incisors and canines



Article# IC6-096
Fig. 67A
Upper Thrid Molars



Article# IC6-097
Fig. 35
Upper Universal



Article# IC6-098
Fig. 35A
Upper Universal



Article# IC6-099
Fig. 79
Lower Thrid Molars



Straight extracting forceps for upper incisors and canines

Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

DEEP ROOT PICK FORCEPS



Article# IC6-100
Fig. 79 T
Lower Third Molars



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines

EXTRACTING FORCEPS



Article# IC6-101
Fig. 44
Apex Forceps Upper Root
Deep Root Pick



Article# IC6-102
Fig. 45
Apex Forceps
Lower Root
Deep Root Pick



Article# IC6-103
Fig. 51
Apex Forceps
Upper Root
Deep Root Pick



Article# IC6-104
Fig. 46XL
Apex Forceps Lower Root
Deep Root Pick



Article# IC6-105
Upper Root
Fig. 44



Article# IC6-106
Lower Root
Fig. 45



Article# IC6-107
Upper Root
Fig. 51



forceps for upper incisors, premolars and roots

■ Инструменты для реставрации
 ■ Instruments pour plombage
 ■ Füllungsinstrumente (mit Kontrawinkel)
 ■ Инструменти для реставрації
 ■ Instrumentos para obturaciones
 ■ ملء الصكوك



Article# IC6-108 F-1 Upper Universal
Article# IC6-109 F-2 Lower Universal
Article# IC6-110 F-3 Upper Anterior
Article# IC6-111 Lower Molar Modify
Article# IC6-112 F-4N Lower Anterior
Article# IC6-123 F-6 Upper 3rd Molar

Article# IC6-115 Extracting Forceps Dual Grip Set of 6 Pcs



Article# IC6-116 Fig. 336 lower premolars deep-gripping
Article# IC6-117 Fig. 1 Upper Anterior
Article# IC6-118 Fig. 35 Upper Universal
Article# IC6-119 Fig. 36 Lower Premolars Deep-Gripping
Article# IC6-120 Fig. 67 Upper Thrid Molars
Article# IC6-121 Fig. 79A Apex Forceps Lower Universal Deep Root Pick

(Deep Gripe Extracting Forceps)



Article# IC6-124 Fig. 321 Lower Universal
Article# IC6-125 Fig. 36 Lower Premolars Deep-Gripping
Article# IC6-126 Fig. 336 lower premolars deep-gripping
Article# IC6-127 Fig. 67 Upper Thrid Molars
Article# IC6-128 Fig. 79 T Lower Thrid Molars



Article# IC6-123 Fig. 321A Lower Root Molars
 • No Fracture • No Tissue Damage • Extraction With Easy

Article# IC6-129 Deep Gripe Extracting Forceps



Article# IC6-114
 Extracting Forceps Dual Grip Set of 6 Pcs



Article# IC6-122
 Extracting Forceps Dual Grip Set of 6 Pcs



Инструменты для реставрации
Instruments pour plombage
Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації
Instrumentos para obturaciones
ملء الصكوك



Fig. 1
upper
laterals
and canines



Article# IC6-130

Set of 10 Pieces

Fig. 1
Fig. 7
Fig. 13
Fig. 17
Fig. 18
Fig. 22
Fig. 51A
Fig. 67
Fig. 74N
Fig. 79

Fig. 7
upper premolars



Fig. 13
lower premolars



Fig. 17
upper molars, right



Fig. 18
upper molars, left



Fig. 22
lower roots



Fig. 51A
upper roots



Fig. 67
upper third molars



Fig. 74N
lower molar



Fig. 79
lower third molars



Extracting forceps for upper roots



Инструменты для реставрации
Instruments pour plombage
Füllungsinstrumente (mit Kontrawinkel)

Інструменти для реставрації
Instrumentos para obturaciones
ملء الصكوك

EXTRACTING FORCEPS



Fig. 1
upper
laterals
and canines



Article# IC6-131
Set Of 10 Pices

Fig. 7
upper premolars



Fig. 13
lower premolars



Fig. 17
upper molars, right



Fig. 18
upper molars, left



Fig. 22
lower roots



Fig. 33A
upper roots



Fig. 51A
upper third molars



Fig. 67A
upper third molars



Fig. 79
lower third molars



Инструменты для реставрации
Instruments pour plombage
Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації
Instrumentos para obturaciones
ملء الصكوك

Pediatric Tooth Extracting Forceps



Article# IC6-132
English Pattern
Klein Fig. 4
Pedodontic Upper root



Article# IC6-133
English Pattern
Klein Fig. 1
Pedodontic, Upper Incisors



Article# IC6-134
English Pattern Klein Fig. 5
Pedodontic Lower Incisors.



Article# IC6-135
English Pattern
Klein Fig. 2
Pedodontic Upper Premolars

Baby Tooth Extraction Forceps




Article# IC6-136
Set of 7 Pieces



Article# IC6-137
English Pattern
Klein Fig. 6
Pedodontic Lower Molars



Article# IC6-138
English Pattern
Klein Fig. 3
Pedodontic Upper Molars Universal



Article# IC6-139
Tooth Forceps Children Set For 7Pcs



Article# IC6-140
English Pattern
Klein Fig. 7
Pedodontic Lower Molars,

- | | | |
|----|---------------|------|
| 01 | Tooth Forceps | Fig. |
| 02 | Tooth Forceps | 1 |
| 03 | Tooth Forceps | 2 |
| 04 | Tooth Forceps | 3 |
| 05 | Tooth Forceps | 4 |
| 06 | Tooth Forceps | 5 |
| 07 | Tooth Forceps | 6 |
| | | 7 |

■ Инструменты для реставрации
 ■ Instruments pour plombage
 ■ Füllungsinstrumente (mit Kontrawinkel)
 ■ Инструменти для реставрації
 ■ Instrumentos para obturaciones
 ■ ملء الصكوك



Article# IC6-141
Witzel Forceps With Gold
Set of 3 Pcs



Article# IC6-145
Witzel Forceps With Hole And Screw Pattern
Set of 3 Pcs



Article# IC6-142
#401(301)
160mm
Lower roots, Universal



Article# IC6-146
#401(301)
160mm
Lower roots, Universal



Article# IC6-143
#400(300)
165mm
Upper roots, Universal



Article# IC6-147
#400(300)
165mm
Upper roots, Universal



Witzel
Article# IC6-144
165mm
Lower roots, Universal



Witze.
Article# IC6-148
165mm
Lower roots, Universal



Article# IC6-149
Fig # 223
(Witzel, 150mm)
Universal for lower roots



Article# IC6-150
Fig # 221
(Witzel, 150mm)
Universal for upper roots



Article# IC6-151
Splinter Forceps
15.5cm



Article# IC6-152
Splinter Forceps
15.5cm 45°



Article# IC6-153
Stieglitz
Root Fragments



Article# IC6-155
Stieglitz Splinter
Forceps
140mm 90°



Article# IC6-156
Stieglitz Splinter Forceps
140mm Straight

Инструменты для реставрации
Instruments pour plombage
Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації
Instrumentos para obturaciones
ملء الصكوك



Article# IC6-157
Extracting Forceps
Set of 10 Pieces
Extracting Forceps Fig. 2
Extracting Forceps Fig. 7
Extracting Forceps Fig. 13
Extracting Forceps Fig. 17
Extracting Forceps Fig. 18
Extracting Forceps Fig. 22
Extracting Forceps Fig. 33
Extracting Forceps Fig. 51A
Extracting Forceps Fig. 67A
Extracting Forceps Fig. 79



Article# IC6-158
Extracting Forceps
Set of 10 Pieces
Extracting Forceps Fig. 2
Extracting Forceps Fig. 7
Extracting Forceps Fig. 13
Extracting Forceps Fig. 17
Extracting Forceps Fig. 18
Extracting Forceps Fig. 22
Extracting Forceps Fig. 33
Extracting Forceps Fig. 51A
Extracting Forceps Fig. 67A
Extracting Forceps Fig. 79



Article# IC6-159
Extracting Forceps (English Pattern)
Set of 13 Pieces
Extracting Forceps Fig. 2
Extracting Forceps Fig. 7
Extracting Forceps Fig. 13
Extracting Forceps Fig. 17
Extracting Forceps Fig. 18
Extracting Forceps Fig. 22
Extracting Forceps Fig. 33
Extracting Forceps Fig. 51A
Extracting Forceps Fig. 67A
Extracting Forceps Fig. 79
Extracting Forceps Fig. 86C
Extracting Forceps Fig. 89
Extracting Forceps Fig. 90



Article# IC6-160
Extracting Forceps (English Pattern)
Set of 10 Pieces
Extracting Forceps Fig. 13S
Extracting Forceps Fig. 22S
Extracting Forceps Fig. 30S
Extracting Forceps Fig. 33S
Extracting Forceps Fig. 37
Extracting Forceps Fig. 38
Extracting Forceps Fig. 39
Extracting Forceps Fig. 39L
Extracting Forceps Fig. 39R
Extracting Forceps Fig. 51S

Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك



Article# IC6-161

Extracting Forceps (English Pattern)

Set of 10 Pieces

Extracting Forceps Fig. 1
Extracting Forceps Fig. 7
Extracting Forceps Fig. 13
Extracting Forceps Fig. 17
Extracting Forceps Fig. 18
Extracting Forceps Fig. 22
Extracting Forceps Fig. 33A
Extracting Forceps Fig. 51A
Extracting Forceps Fig. 67A
Extracting Forceps Fig. 79



Article# IC6-162

Extracting Forceps (Anatomical Handle)

Set of 10 Pieces

Extracting Forceps Fig. 1
Extracting Forceps Fig. 7
Extracting Forceps Fig. 13
Extracting Forceps Fig. 17
Extracting Forceps Fig. 18
Extracting Forceps Fig. 22
Extracting Forceps Fig. 33A
Extracting Forceps Fig. 51A
Extracting Forceps Fig. 67A
Extracting Forceps Fig. 79



Article# IC6-163

Extracting Forceps (Hole Pattern Handle)

Set of 10 Pieces

Extracting Forceps Fig. 1
Extracting Forceps Fig. 7
Extracting Forceps Fig. 13
Extracting Forceps Fig. 17
Extracting Forceps Fig. 18
Extracting Forceps Fig. 22
Extracting Forceps Fig. 33A
Extracting Forceps Fig. 51A
Extracting Forceps Fig. 67A
Extracting Forceps Fig. 79

 Инструменты для реставрации

🇫🇷 Instruments pour plombage

- Füllungsinstrumente (mit Kontrawinkel)

Інструменти для реставрації

 Instrumentos para obturaciones

ملء الصكوك



Article# IC6-164

Extracting Forceps English adult Adult Set of 10 Pcs

01	Tooth Forceps	Fig. 1
02	Tooth Forceps	7
03	Tooth Forceps	13
04	Tooth Forceps	17
05	Tooth Forceps	18
06	Tooth Forceps	22
07	Tooth Forceps	33
08	Tooth Forceps	51
09	Tooth Forceps	67
10	Tooth Forceps	79



Article# IC6-165

Tooth Forceps Set of 8 Pcs

01 Tooth Forceps
02 Tooth Forceps
03 Tooth Forceps
04 Tooth Forceps
05 Tooth Forceps
06 Tooth Forceps
07 Tooth Forceps
08 Tooth Forceps



Article# IC6-166

Extracting Forceps English adult Adult Set of 10 Pcs

01	Tooth Forceps	Fig. 1
02	Tooth Forceps	7
03	Tooth Forceps	13
04	Tooth Forceps	17
05	Tooth Forceps	18
06	Tooth Forceps	22
07	Tooth Forceps	33
08	Tooth Forceps	51
09	Tooth Forceps	67
10	Tooth Forceps	79



Article# IC6-167
Set of 37 Pcs

- 01 - Extracting Forces Fig.1
- 02 - Extracting Forces Fig.7
- 03 - Extracting Forces Fig.13
- 04 - Extracting Forces Fig.17
- 05 - Extracting Forces Fig.18
- 06 - Extracting Forces Fig.18A
- 07 - Extracting Forces Fig.19
- 08 - Extracting Forces Fig.20
- 09 - Extracting Forces Fig.22
- 10 - Extracting Forces Fig.29
- 11 - Extracting Forces Fig.33
- 12 - Extracting Forces Fig.51
- 13 - Extracting Forces Fig.51A
- 14 - Extracting Forces Fig.67A
- 15 - Extracting Forces Fig.79
- 16 - Extracting Forces Fig.136
- 17 - Extracting Forces Fig.137
- 18 - Extracting Forces Fig.53 L
- 19 - Extracting Forces Fig.53 R
- 20 - Extracting Forces Fig.65
- 21 - Root Elevators Fig.1B
- 22 - Root Elevators Fig.3B
- 23 - Root Elevators Fig.5B
- 24 - Root Elevators Fig.4B SA Apexo
- 25 - Root Elevators Fig.5B SA Apexo
- 26 - Cotton Tweezers
- 27 - Mouth Mirrors With Handles Fig.5
- 28 - Explorers



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines



Article#

IC6-168 Available in Gold Coated Set of 25 Pcs

IC6-169 Available in Polish Set of 25 Pcs

- 01 - Extracting Forceps Fig.13
- 02 - Extracting Forceps Fig.17
- 03 - Extracting Forceps Fig.18
- 04 - Extracting Forceps Fig.19
- 05 - Extracting Forceps Fig.20
- 06 - Extracting Forceps Fig.22
- 07 - Extracting Forceps Fig.29
- 08 - Extracting Forceps Fig.33
- 09 - Extracting Forceps Fig.51
- 10 - Extracting Forceps Fig.51A
- 11 - Extracting Forceps Fig.67A
- 12 - Extracting Forceps Fig.79
- 13 - Root Elevators Fig.1B
- 14 - Root Elevators Fig.3B
- 15 - Root Elevators Fig.4B SA Apexo
- 16 - Root Elevators Fig.5B SA Apexo
- 17 - Cotton Tweezers
- 18 - Mouth Mirrors With Handles
- 19 - Explorers
- 20 - Cement Spatulas
- 21 - Filling Instruments
- 22 - Filling Instruments
- 23 - Filling Instruments
- 24 - Excavators
- 25 - Dental Aspirated Syringes



Article# IC6-170
Luxating Elevator
Straight 3mm



Article# IC6-171
Luxating Elevator
Curved 3mm



Article# IC6-183
Goldman Fox
Scissor Cvd



Article# IC6-184
Crile Wood
Needle Holder TC



Article# IC6-185
Atraumatic Extraction Kit
of 16 Pcs

Atraumatic Extraction Kit



Article# IC6-172 1 Periotome



Article# IC6-173 2 Periotome



Article# IC6-174 9 Molt Periosteal Black



Article# IC6-175 Vista Tunneling Elevator 1



Article# IC6-176 Scalpel Handle 5 Straight



Article# IC6-177 1 Atraumatic Extraction Forcep



Article# IC6-178 35N Atraumatic Extraction Forcep



Article# IC6-179 36N Atraumatic Extraction Forcep



Article# IC6-180 Cawood Minnesota Cheek Retractor 16cm



Article# IC6-181 Urban Cheek Retractor 14cm



Article# IC6-182 Adson Brown Tissue Plier

Щипцы экстракционные (американская форма)
Daviers - forme américaine
Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)
Pinzeper per estrazione - forma americana
ستخراج الملقط (النمط الأمريكي)



Article# IC6-186
1STD
Winter
upper centrals and canines

Article# IC6-187
Fig. 1
upper incisors
and roots

Article# IC6-188
Cryer
Fig. 150

Щипцы экстракционные (американская форма)
 Daviers - forme américaine
 Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)
 Pinzeper per estrazione - forma americana
 استخراج الملقط (النمط الأمريكي)



Straight extracting forceps for upper incisors and canines




Article# IC6-190
 Fig. 62
 incisors, premolars,
 deciduous teeth



Article# IC6-191
 # 10S
 upper roots



Article# IC6-192
 # 17
 lower molars



Article# IC6-193
 Harris
 Fig. 18R



Article# IC6-194
 Harris
 Fig. 18L

Article# IC6-189
 Fig. 150A
 Cryer
 upper incisors,
 premolars, roots



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines

Щипцы экстракционные (американская форма)
Daviers - forme américaine
Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)
Pinzeper per estrazione - forma americana
ستخراج الملقط (النمط الأمريكي)



Article#
IC6-195
Fig. 210
upper third molars



Article#
IC6-196
Fig. 53R
upper molars, right



Article#
IC6-197
Fig. 53L
upper molars, left



Article#
IC6-198
NEVIUS
Fig. 88R



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines

Щипцы экстракционные (американская форма)
 Daviers - forme américaine
 Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)
 Pinzeper per estrazione - forma americana
 استخراج الملقط (النمط الأمريكي)



Article# IC6-199
 Fig. 24
 upper molars



Article# IC6-200
 # 210S
 upper third molars



Article# IC6-201
 # 210H
 upper third molars



Article# IC6-202
 # 286
 Upper Roots,
 Premolars, Incisors



Article# IC6-203
 # 99C
 Kells
 upper incisors,
 canines, premolars



Article# IC6-204
 # 222
 lower third molars



Article# IC6-205
 # 286
 Lower Molar

Щипцы экстракционные (американская форма)
Daviers - forme américaine
Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)
Pinzeper per estrazione - forma americana
ستخراج الملقط (النمط الأمريكي)



Article# IC6-206
Nevius
Fig. 88L
upper molars, left



Article# IC6-207
Fig. 65
upper incisors
and roots



Article# IC6-208
Parmly Fig. 32A
upper premolars
and roots




Article# IC6-209
Parmly Fig. 32
upper premolars
and molars

Щипцы экстракционные (американская форма)
 Daviers - forme américaine
 Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)
 Pinzeper per estrazione - forma americana
 استخراج الملقط (النمط الأمريكي)

EXTRACTING FORCEPS



Article# IC6-210
 COHEN
 universal

Article# IC6-211
 Fig. 203
 lower incisors
 and roots

Article# IC6-212
 Tomes Fig. 69
 upper and lower
 roots

Щипцы экстракционные (американская форма)
Daviers - forme américaine
Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)
Pinzeper per estrazione - forma americana
ستخراج الملقط (النمط الأمريكي)



Article# IC6-213
Fig. 23
lower molars

Article# IC6-214
Fig. 16
lower molars

Article# IC6-215
Physic Fig. 5
lower third molars



Article# IC6-216
Harris Fig. 15
lower molars



Article# IC6-217
Fig. 103

Щипцы экстракционные (американская форма)
 Daviers - forme américaine
 Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)
 Pinzeper per estrazione - forma americana
 استخراج الملقط (النمط الأمريكي)

EXTRACTING FORCEPS



Article# IC6-218
 Hull
 Fig. 101
 Deciduous teeth



Article# IC6-219
 Cryer
 Fig. 150S
 upper teeth,
 for children



Article# IC6-220
 Fig. 150 AS

Щипцы экстракционные (американская форма)
Daviers - forme américaine
Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)
Pinzeper per estrazione - forma americana
ستخراج الملقط (النمط الأمريكي)



Article# IC6-221
Cryer
Fig. 151
lower teeth,
for children

Article# IC6-222
Fig. 151 S
lower incisors,
premolars, roots

Article# IC6-223
Fig. 151 AS

Щипцы экстракционные (американская форма)
 Daviers - forme américaine
 Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)
 Pinzeper per estrazione - forma americana
 ستخراج الملقط (النمط الأمريكي)



Article# IC6-224
16S
 lower molars for children



Article# IC6-225
17S
 lower molars



Article# IC6-226
23 S
 Lower molars



Article# IC6-227
27
 lower molars



Article# IC6-228
18LS
 upper molars left



Article# IC6-229
18RS
 upper molars right



Article# IC6-230
85 Tomas
 lower canines, premolars and molars



Article# IC6-231
85A Tomas
 lower canines, premolars and molars



Article# IC6-232
69S
 upper and lower roots



Article# IC6-233
99A Kells
 upper incisors, canines, premolars



Article# IC6-234
217
 lower molars

Щипцы экстракционные (американская форма)
Daviers - forme américaine
Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)
Pinzeper per estrazione - forma americana
ستخراج الملقط (النمط الأمريكي)



Article# IC6-235
#10 SK
Upper Molars Sup



Article# IC6-236
#17 SK Universal
Lower Molars



Article# IC6-237
#23 SK Universal
Upper Molars



Article# IC6-238
69 SK Tomas
upper and lower roots



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines



Article# IC6-239
99 SK
upper incisors,
canines, premolars



Article# IC6-240
150 SK
upper incisors,
premolars, roots



Article# IC6-241
151 SK
lower incisors,
premolars, roots



Article# IC6-242
222 SK
lower third molars

- Щипцы экстракционные (американская форма)
- Давиеры - forme americane
- Férceps de extracción - forma americana
- Щипці екстракційні (американська форма)
- Pinzeper per estrazione - forma americana
- ستخراج الملقط (النمط الأمريكي)



Article# IC6-243
Fig. 1
Upper Anterior



Article# IC6-244
Lower Molar
Fig. 23 Cow Horn



Article# IC6-245
Fig. 88L
Upper Molar Left



Article# IC6-246
Fig. 88R
Upper Molar Right



Article# IC6-247
Fig. 150
Upper Universal



Article# IC6-248
Fig. 150S
Upper Universal



Article# IC6-249
Fig. 151
Lower Universal



Article# IC6-250
Fig. 151S
Lower Universal



Article# IC6-251
Fig. 106
Lower Incisors



Article# IC6-252
Fig. 174
Extra Narrow Lower Anterior
Apical Retention Forceps

Шприцы стоматологические карпульные
Seringues à usage dentaire
Dentalspritzen

Шприці стоматологічні карпульні
Jeringas dentales
المحاقن الأسنان



Type C
Article#
IC6-253 1.8ml
IC6-254 2.2ml



Type C
Article#
IC6-255 1.8ml
IC6-256 2.2ml



Type C
Article#
IC6-257 1.8ml
IC6-258 2.2ml




Article# IC6-259
Aspirated Dental Cartridge
Syringe 1.8 ml




Article# IC6-260
Aspirated Dental Cartridge
Syringe Folding Style 1.8 ml



Article# IC6-261
Aspirated Dental Cartridge
Syringe Folding Style 1.8 ml



Euro Breech loading
Article#
IC6-262 1.8ml
IC6-263 2.2ml



EuroBreech loading
Article#
IC6-264 1.8ml
IC6-265 2.2ml



Article#
IC6-266 1.8ml
IC6-267 2.2ml



Article# IC6-268
Aspirated Dental Cartridge
Syringe Folding Style 1.8 ml



Aspirating Syringe
Article#
IC6-269 1.8mm
IC6-270 2.2mm



Article#
IC6-271 1.8mm
IC6-272 2.2mm



Aspirating Syringe
Article#
IC6-273 1.8mm
IC6-274 2.2mm



2-Ring
Article#
IC6-275 1.8ml
IC6-276 2.2ml



Aspirating Syringe
Article#
IC6-277 1.8mm
IC6-278 2.2mm



2-Ring
Article#
IC6-279 1.8ml
IC6-280 2.2ml



Cartridge Syringe
Article#
IC6-281 1.8mm
IC6-282 2.2mm



Cartridge Syringe
Article#
IC6-283 1.8mm
IC6-284 2.2mm



Cartridge Syringe
Article#
IC6-285 1.8mm
IC6-286 2.2mm



Article#
IC6-287 1.8ml
IC6-288 2.2ml



Manual Aspirating Syringe
Article#
IC6-289 1.8ml
IC6-290 2.2ml



Manual Aspirating Syringe
Article#
IC6-291 1.8ml
IC6-292 2.2ml



Quatro
Article#
IC6-293 1.8ml
IC6-294 2.2ml



Trio
Article#
IC6-295 1.8ml
IC6-296 2.2ml



Article#
IC6-297 1.8ml
IC6-298 2.2ml



Article#
IC6-299 1.8ml
IC6-300 2.2ml

Шприцы стоматологические карпульные
Seringues à usage dentaire
Dentalspritzen

Шприці стоматологічні карпульні
Jeringas dentales
المحاقن الأسنان

EXTRACTING FORCEPS



Disc System
Article#
IC6-301 1.8ml
IC6-302 2.2ml



Disc System
Article#
IC6-303 1.8ml
IC6-304 2.2ml



Disc System
Article#
IC6-305 1.8ml
IC6-306 2.2ml



Projection System
Article#
IC6-307 1.8ml
IC6-308 2.2ml



Projection System
Article#
IC6-309 1.8ml
IC6-310 2.2ml



Article# IC6-311
Dental Anesthesia Syringe
1.8 ml (in Multi Colors)



Article# IC6-312
Dental Anesthesia Syringe
1.8 ml (in Multi Colors)



Article# IC6-313
Dental Anesthesia Syringe
1.8 ml (in Multi Colors)



Article# IC6-314
Suitable For Use With All
1.8ml Dental Cartridges




Article# IC6-315
Suitable For Use With All
1.8ml Dental Cartridges



Article#
IC6-316 1.8ml
IC6-317 2.2ml
Angled



Ject Ultima
Article#
IC6-318 1.8ml
IC6-319 2.2ml
Injects 0.05 ml with every click Supplied with
angled and straight tubes



Ject angled
Article#
IC6-320 1.8ml
IC6-321 2.2ml



Ject straight
Article#
IC6-322 1.8ml
IC6-323 2.2ml



Swedish Type
Article#
IC6-324 1.8ml
IC6-325 2.2ml



Article# IC6-326
1001
Intraligamental
Anesthesia
Syringe in
Black Colour



Article#
IC6-327 1.8ml
IC6-328 2.2ml



Aluminium - Dark
Article#
IC6-329 1.8 ml
IC6-330 2.2 ml



Aluminium - Silver
Article#
IC6-331 1.8 ml
IC6-332 2.2 ml



Aluminium - Orange
Article#
IC6-333 1.8 ml
IC6-334 2.2 ml



Aluminium - Green
Article#
IC6-335 1.8 ml
IC6-336 2.2 ml



Aluminium - Blue
Article#
IC6-337 1.8 ml
IC6-338 2.2 ml